

CNW Solvent Filter

CNW Solvent Filter is designed to solvent rapid filtration and degassing, suitable for filtering and degassing of the HPLC mobile phase solvent, Can prolong the service life of instrument and the chromatographic column, improve the detection accuracy; Used in gravimetric analysis, trace analysis, trace analysis, colloid separation and sterile in the laboratory.



Compare to traditional Solvent Filter , the new design solve following questions:

1. Larger filter cup capacity, avoid multiple adding;
2. Removes clamp, avoid potential leaking;
3. Substitute Frosted Seal to Thread fastening, avoid fusion after long time usage.
4. Substitute Quartz sand core to PTFE filter plate, avoid hard cleaning.
5. Substitute Conical flask to GL45 ISO Bottles, avoid solvent transfer after filter.



Cat.No.	Description	Packing
2.PSF010.0001	1000mL Vacuum Filter Holder (NOT include GL45 ISO Bottle)	1 set